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INDICATIONS

WINN'SAVE® devices for restraining or securing patients in a bed are subject to medical prescription. These products are designed and used for safety reasons for beds in the following cases:

- Physical restraint: to restrict the movements of an agitated and/or aggressive patient presenting cognitive or mental problems, to stop him/her falling and to protect staff or the patient from him/herself.
- Securing: protecting a confused patient from getting up unauthorised for medical reasons and for administering treatment.

PRECAUTIONS OF USE

This product must be attached by qualified personnel according to the manufacturer's instructions by following the safety regulations and procedures in line with the best practice recommendations applicable in health centres and regarding home care.

- Follow the instructions for use: restraining or securing in a chair depending on the indication, choice of size depending on waist circumference, adjustment of the length of the removable pelvic strap, comfortable tightness that does not restrict normal breathing at rest, specific surveillance of the skin in direct contact (rubbing, maceration).
- Before use, check that the characteristics of the restraining device are intact (eyelets, loops, seams, buckles, pelvic strap).
- Always attach the waist restraint to the pelvic strap using the side straps, unless contraindicated by the prescriber.
- In the event of direct contact with a patient in isolation for a skin condition (contact precautions), the device will be considered for use by said patient only, unless reuse is authorised after treatment by the hospital hygiene department.
- Use only the closing systems provided or approved by the manufacturer.

SPECIFIC SURVEILLANCE

- Prohibit any alteration or repair to the restraining or securing devices.
- Check that an agitated patient cannot access or handle the closing systems.
- Check the restraining devices used and stocked for the intention of physical restraint at regular intervals. They must be clean, complete, in good condition and kept such that they can be used swiftly.
- Keep unused devices out of patients' reach.
- Check that cleaning and maintenance are carried out in accordance with the instructions given in the user's manual.
- Check that the buckle for adjusting the length of the removable pelvic strap is protected between the 2 rear fabrics of the waist restraint.

Reference	Designation	Measurements	Closing system	
			Magnet	Buckle and Velcro
VCBTLT002-10	Waist and pelvic restraint for the bed with side straps (S)	45-85 cm*	X	
VCBTLT002-20	Waist and pelvic restraint for the bed with side straps (M)	70-118 cm*	X	
VCBTLT002-21	Waist and pelvic restraint for the bed without side straps (M)	70-118 cm*	X	
VCBTLT002-30	Waist and pelvic restraint for the bed with side straps (L)	94-158 cm*	X	
VCBTLT002-31	Waist and pelvic restraint for the bed without side straps (L)	94-158 cm*	X	
VCBTLT003-01	Optional: bed rail for ankle trap	-	X	
VCBTLT004-01	Chest harness extention	-	X	
VCBTLT004-20	Chest harness restraint	20-50 cm	X	
VCBTLT006-00	Universal strap restraints for wrists and ankles (x2)	5-40 cm**	X	
VCBTLT011-00	One-size strap restraints for wrists (x2)	10-30 cm**	X	
VCNTLT007-00	Universal strap with buckle (x2)	5-40 cm**		X
VCNTLT010-20	Waste strap (Buckles, Velcro) - (M)	65-85 cm*		X
VCNTLT010-30	Waste strap (Buckles, Velcro) - (L)	85-120 cm*		X

*Measurement corresponding to waist circumference

**Measurement corresponding to ankle and wrist circumference

Reference	Designation
61-0605-15	Complete magnetic closing with magnetic key
61-0605-30	Complete magnetic closing with no magnetic key
61-0605	Magnetic key for magnetic closing

INSTRUCTIONS

Magnetic closing

Put the stainless steel tip through the eyelet hole in the fabric strap. Insert the strap eyelet intended to close the restraining device. Close by clicking the button on the tip. The magnetic key will open the device by removing the button.


Precautions: when installing a restraining device, check that the closing system is correctly locked. To facilitate opening, release the pressure applied to the button.

Device for restraining or securing the patient's chest to the bed


Precautions: the waist restraint straps are fastened in the middle of the bed's long sides. Tighten until crushing the edge of the mattress. Lock the moving functions of the bed if possible (back and leg rests). Adjust the length of the removable pelvic strap before use according to the patient's morphology.

Chest restraint with or without extensions

Precautions: In the case of a medical bed with adjustable height, fasten the device to the mobile parts that move with the patient. In the case of a non-medical bed, fasten the device to the bed base or frame. Lock the moving functions of the bed if possible (back and leg rests).


Device for restraining or securing wrists

Secure wrists on the waist restraint strap or to the fixed side part of the bed.

Secure wrists using the buckle and fasten to the fixed structure of the bed.

Ankle restraints

Restrain ankles by fastening to the fixed structure of the bed rail. In the case of severe agitation, the bed rail can be secured using an extension secured to the middle of the rail and bed frame. Lock the moving functions of the bed base.

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